






















Shortcut Instructions for Desktop Scanners





1, Model operation

 Rotate Model	 Scale Model	 Translate Model	 Add Scan
--	---	---	--



2, The shortcut in scan process

 Start Scan	Space	 Articulator Dynamic Scan	P
 Finish	Space	 High Quality	D
 Cancel Scan	Esc	HDR HDR Mode	H
 AI Scan	A	 Choose Detailed Area	H
 Texture Scan	C	 Reduce High Brightness	R
 Show Texture	C	 Abutment Scan	B
 Multi Path Scan/Quick bite scan	P	 Unsectioned Die	U
 Swap Jaws	W		
 Import Model	I		
 Edit	E		
 New Order	N		

3, The shortcut in pre-design process

 Edit	E
 Fill Holes	F
 Show Texture	C
 Smooth	M

4, The shortcut for entire software

 Screenshot	S
 Shortcut Instructions	K