

SHINING 3D DENTAL

**3rd Party Material Packages available for AccuFab
Printers(D1s,L4D)**

Brand	Material	Layer Thickness	AccuFab-D1s	AccuFab-L4D
Keystone	KeyCast	0.05	✓	✓
	KeyDenture Try-In	0.05	✓	✓
	KeyGuide	0.1	✓	NA
	KeyMask	0.1	✓	✓
	KeyModel	0.1	✓	✓
	KeyOrtho IBT	0.1	✓	✓
	KeyOrthoModel	0.1	✓	✓
	KeySplint Hard	0.1	✓	✓
	KeySplint Soft	0.1	✓	✓
	KeySplint Soft fast (Only for V5 fast resin tank)	0.1	✓	NA
	KeyTray	0.1	✓	
pro3dure	GR-10 guide	0.1	✓	NA
	GR-11 tray	0.1	✓	
	GR-12 cast	0.05	✓	
	GR-14.1 denture	0.075	✓	✓
	GR-14.2 deep pink	0.075	✓	✓
	GR-14.2 light pink	0.075	✓	NA
	GR-14.2 orange pink	0.075	✓	✓
	GR-16 Xray	0.05	✓	✓
	GR-17 temporary	0.05	✓	NA
	GR-17 temporary A1	0.075	NA	✓
	GR-17.1 temporary lt	0.05	✓	NA
	GR-18.1 IB	0.1	✓	✓
	GR-19 OA	0.1	✓	✓
	GR-21 Try-In	0.05	✓	✓

	GR-22 flex	0.1	✓	✓
	HI GR-14.2 orange pink fast (Only for V5 fast resin tank)	0.075	✓	NA
HARZ Labs	DENTAL BLEACH	0.05	✓	✓
	DENTAL PINK	0.075	✓	✓
	DENTAL SAND A1A2	0.05	✓	✓
DETAX	FREEPRINT cast 2.0 385	0.05	✓	✓
	FREEPRINT denture 385 translucent pink	0.05	✓	✓
	FREEPRINT gingiva 385	0.05	✓	✓
	FREEPRINT IBT 385	0.1	✓	✓
	FREEPRINT model T 385	0.1	✓	✓
	FREEPRINT ortho 385	0.075	✓	✓
	FREEPRINT ortho 405	0.05	NA	✓
		0.075	✓	NA
	FREEPRINT splint 2.0 385	0.1	✓	✓
	FREEPRINT temp 385 A2	0.05	✓	✓
dentona	optiprint clara transparent	0.1	✓	✓
	optiprint laviva light pink	0.1	✓	✓
	optiprint lumina A2	0.05	✓	✓
	optiprint zero red	0.05	✓	✓
	optiprint gingiva 385	0.1	NA	✓
	optiprint guide 405	0.075		✓
	optiprint IBT 385	0.1		✓

	optiprint match	0.05		✓
	optiprint model align	0.1		✓
	optiprint model golden brown	0.075		✓
	optiprint model grey	0.075		✓
	optiprint model ivory	0.075		✓
	optiprint tray HR	0.075		✓
Liqcreate	Clear Impact	0.05	NA	✓
		0.1		✓
	Composite-X	0.05		✓
		0.1		✓
	Deep Blue	0.05		✓
		0.1		✓
	Dental Model Pro Beige	0.05		✓
		0.1		✓
	Dental Model Pro Grey	0.05		✓
		0.1		✓
	Flexible-X	0.05		✓
		0.1		✓
	Gingiva Mask	0.05		✓
		0.1		✓
	Premium Black	0.05		✓
		0.1		✓
	Premium Flex	0.1		✓
	Premium Model	0.05		✓
		0.1		✓
	Premium Tough	0.1		✓
	Premium White	0.05		✓
		0.1		✓
	Stone Coal Black	0.05		✓
		0.1		✓

	Strong-X	0.05		✓
		0.1		✓
	Tough-X	0.05		✓
		0.1		✓
	Wax Castable	0.025		✓
		0.05		✓
		0.1		✓
Gorky Liquid	Castable	0.05	✓	✓
	Dental Base	0.05	✓	✓
	Dental Castable	0.05	✓	✓
	Dental Crown A1A2	0.05	✓	✓
	Dental Model	0.05	✓	✓
	Dental Surgical	0.05	✓	✓
	Flex	0.1	✓	✓
	Force	0.1	✓	✓
	Reactive	0.1	✓	✓
	Simple	0.1	✓	✓
Dreve	BioTec	0.1	✓	✓
	FotoDent denture	0.05	✓	✓
	FotoDent IBT	0.1	✓	✓
	FotoDent tray	0.1	✓	✓
BASF	Ultracur 3D DM 2304	0.05	NA	✓
		0.1		✓
	Ultracur 3D DM 2505	0.05	✓	NA
	Ultracur 3D EL 4000	0.05	NA	✓
		0.1		✓
	Ultracur 3D EL 4000 B	0.05		✓
		0.1		✓
	Ultracur 3D EL150	0.05	✓	✓
		0.1	NA	✓
	Ultracur 3D FL60	0.05		✓

		0.075	✓	✓
		0.1	NA	✓
	Ultracur 3D RG 1100	0.1		✓
	Ultracur 3D RG 35	0.05	✓	✓
		0.075		NA
		0.1		✓
	Ultracur 3D RG 35 B	0.05	NA	✓
		0.1		✓
	Ultracur 3D RG 50	0.05	✓	NA
		0.1	NA	✓
	Ultracur 3D ST 7500 G	0.05		✓
		0.1		✓
	Ultracur 3D ST 80	0.05	✓	✓
		0.1	NA	✓
	Ultracur 3D ST 80 B	0.05		✓
		0.1		✓
	Ultracur 3D ST 80 G	0.05		✓
	Ultracur 3D ST 80 W	0.05		✓
		0.		✓
	Ultracur 3D ST 45 B	0.05	✓	✓
		0.075	✓	✓
		0.1	✓	✓
	Ultracur 3D ST 45 M	0.05	✓	NA
		0.075	✓	
		0.1	✓	
DMG	LuxaPrint Cast	0.05	✓	✓
	LuxaPrint Model GRY	0.05	✓	✓
		0.1	✓	✓
	LuxaPrint Model IVY TRA	0.1	✓	✓
	LuxaPrint Ortho	0.05	✓	✓

	LuxaPrint Ortho Plus	0.05	✓	✓
	LuxaPrint Tray	0.1	✓	✓
NextDent	C&B MFH N1	0.05	✓	✓
	C&B MFH N1 fast (Only for V5 resin tank)	0.1	✓	NA
	C&B MFH N1.5	0.05	✓	✓
	C&B MFH N2	0.05	✓	✓
	C&B MFH N2.5	0.05	✓	✓
	Denture 3D+ Dark Pink	0.05	✓	✓
	Denture 3D+ Light Pink	0.05	✓	✓
	Denture 3D+ Opaque Pink	0.05	✓	✓
	Denture 3D+ Red Pink	0.05	✓	✓
	MFH-N1.5	0.05	NA	✓
	MFH-N2.5	0.05		✓
	Ortho Flex	0.1	✓	✓
Graphy	BR-23-5	0.05	✓	✓
		0.1	✓	✓
	EP-400	0.05	✓	✓
		0.1	✓	✓
	S-100M Beige	0.05	✓	✓
		0.1	✓	✓
	S-100M Gray	0.05	✓	✓
		0.1	✓	✓
	SC-130	0.05	✓	✓
		0.1	✓	✓
	SG-100	0.05	✓	✓
		0.1	✓	✓

	TC-80DP	0.05	✓	✓
		0.1	✓	✓
	TC-85DAC	0.05	✓	✓
		0.1	✓	✓
	TFDH	0.05	✓	✓
		0.1	✓	✓
	THD	0.05	✓	✓
		0.1	✓	✓
SeramicCAM	Gingiva Mask	0.05	NA	✓
		0.1		✓
	ProModel	0.05		✓
		0.075		✓
		0.1		✓
VOCO	V-Print dentbase (Otoflash)	0.1	NA	✓
	V-Print SG (Otoflash)	0.05		✓
		0.1		✓
	V-Print Splint (Otoflash)	0.1		✓
	V-Print Splint Comfort (Otoflash)	0.1		✓
Bego	VarseoSmile Crown plus A1	0.05	NA	✓
		0.1		✓
	VarseoSmile Crown plus A1 fast (Only for V5 fast resin tank)	0.05	✓	NA
		0.1	✓	
	VarseoSmileTemp A2	0.05	NA	✓

NA = not available